


Alarm Signal Status Unit

(Plant To Alarm Interface)

for the

The SDX-15 Alarm System



Issue 7 12-02-2023

Alarm Signal Status Unit

About this manual.



When you see this symbol, the associated text in **bold type** refers to something which may cause danger or damage.

The Alarm Signal Status Unit is a low cost means of complying with C11 & HTM02 requirements for indication of alarm status in plant rooms, manifold rooms & VIE compounds. It is mounted in a clear fronted polycarbonate enclosure, protected to IP65, suitable for use in external locations without further protection. It consists of the following indication:-

- Alarm conditions 4no
- Normal
- Power On
- System fault

The ASSU also incorporates:

- A one gas transmitter suitable for use with the SDX-15 alarm system
- Four relays, providing volt-free normally closed contacts for connection to other site systems
- A battery reserve to power the unit for a minimum of 4 hours in the event of power loss.

The ASSU monitors the wiring between the plant (or manifold) and itself via TB4 termination boards mounted in or near the contact source, checking for short or open circuits of the cable. If a fault is detected on the contact lines, an alarm condition is transmitted for the condition relating to the faulty line, the system fault lamp will flash and the 5mm round Yellow LED (on Power Supply PCB) will flash.

The system fault lamp will also flash if the ASSU fails to receive signals from a central alarm, here the 5mm round Red LED (on Power Supply PCB) will flash.

In the event of a mains power failure

the system fault lamp will also flash and the 5mm round Green LED (on Power Supply PCB) will go off..

Selecting channel F inhibits the system fault resulting from a loss of signals from a central alarm.

The channel on which the service is to be transmitted is selected with a 16 way rotary switch. Each service on the system is allocated a channel when the system is initially set up, this being entered on the log sheet. The alarm contacts on the plant (or manifold) are connected via the termination board to the input terminals on the ASSU as follows:-

C	Common
1	First condition
2	Second condition
3	Third condition
4	Pressure fault

Any condition not transmitted from this AS SU must be terminated with a 56k resistor to set the condition to normal if the condition is not to be used, or with a 1k8 resistor if the condition is to be transmitted from another location. The resistors are used to prevent a system fault due to short or open circuit.

Note that if a resistor is fitted, the condition must NOT be connected to the termination board. Resistor codes:-

1k8 brown grey red silver, gold or red
56k green blue orange silver, gold red




It is important to ensure that any condition is only transmitted from one location although other conditions on the service may be transmitted from other transmitters.



WARNING. This equipment is not suitable for connection to an IT power system. A readily accessible means of disconnecting the supply must be provided. The maximum prospective fault current must not exceed 1500 amps.

A 240 volt, 50/60 Hertz supply is required, which is connected to terminals under the mains terminal cover. The 2 core inter-panel wiring is connected to signal terminals 1 & 2. The cable screen must be connected to the "S" terminal.

NUMBER OF GASES TO BE TRANSMITTED (ALWAYS SET TO 1)	THE ALARM SIGNAL STATUS UNIT CAN TRANSMIT A SYSTEM FAULT SIGNAL TO THE ALARM SYSTEM.THESE SETTINGS DETERMINE WHERE THE SIGNAL WILL BE TRANSMITTED.																
1 2 3 4	<table border="1"> <thead> <tr> <th>NORMAL</th> <th>WINDOW</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1 2 </td> </tr> <tr> <td>2</td> <td>2 </td> </tr> <tr> <td>3</td> <td>3 </td> </tr> <tr> <td>4</td> <td>4 </td> </tr> </tbody> </table>	NORMAL	WINDOW	1	1 2	2	2	3	3	4	4	<table border="1"> <thead> <tr> <th>CHANNEL</th> </tr> </thead> <tbody> <tr> <td> NO SYSTEM FAULT TRANSMISSION</td> </tr> <tr> <td> CHANNEL D</td> </tr> <tr> <td> CHANNEL E</td> </tr> <tr> <td> NO SYSTEM FAULT TRANSMISSION.USE THIS SETTING WHEN TERMINATION BOARDS ARE NOT USED.</td> </tr> </tbody> </table>	CHANNEL	NO SYSTEM FAULT TRANSMISSION	CHANNEL D	CHANNEL E	NO SYSTEM FAULT TRANSMISSION.USE THIS SETTING WHEN TERMINATION BOARDS ARE NOT USED.
NORMAL	WINDOW																
1	1 2																
2	2																
3	3																
4	4																
CHANNEL																	
NO SYSTEM FAULT TRANSMISSION																	
CHANNEL D																	
CHANNEL E																	
NO SYSTEM FAULT TRANSMISSION.USE THIS SETTING WHEN TERMINATION BOARDS ARE NOT USED.																	



DECLARATION OF CONFORMITY

EU EMC Directive 2014/30/EU
UK Electromagnetic Compatibility
Regulations 2016 SI 2016 No.1091
The Low Voltage Directive 2006/95/EG
Medical Devices Directive 93/42/EEC

Shire Controls Ltd
Studio 3, Channoeks Farm
Gilston, Harlow
Essex, CM20 2RL

Tel. 01279 434399

shirecontrols.com
shiremedical.co.uk

Manufacturer

Shire Controls Ltd
Studio 3 Channoeks Farm,
Gilston, Harlow
Essex CM20 2RL
United Kingdom

Product Type

Part Of SDX-15 Plant Alarm System
Which Includes Plant To Alarm Interface

Year of manufacture 2023

Standards used (2019)

To which this declaration relates is in conformity
with the following standards:

EN60601-1-2 4th Edition 2015

Emissions Standard for Medical Equipment

EN55011, Class A 2016 + A1:2017

Emissions Standard for ISM Equipment

EN60601-1-2 4th Edition 2015

Immunity Standard for Medical Equipment

EN61000-4-2 2009

ESD Requirements

EN61000-4-3 2006 + A1 + A2

Radiated Susceptibility

EN61000-4-4 2016

Electrical Fast Transient Burst Requirement

EN61000-4-5 2017

Surges Requirements

EN61000-4-6 2014

Conducted Susceptibility

EN61000-4-8 2012

Magnetic Field Immunity

EN61000-4-11 2017

Voltage Dips and Interruptions



This work is licensed under the Creative Commons Attribution-NonCommercial 2.5 License. To view a copy of this license, visit <http://creativecommons.org/licenses/by-nc/2.5/> or send a letter to Creative Commons, 543 Howard Street, 5th Floor, San Francisco, California, 94105, USA.

Authorised representative

Director

Signature

**UK
CA**